

# Info Sheet: Glycan/Glycoprotein analysis

## Identification of glycoproteins

### General description:

The N-glycosylation occurs at asparagine residues in the consensus sequence N-X-S/T (X≠P). The workflow consists of:

- Glycoprotein enrichment with specific protocols
- Protein digestion (Trypsin) to generate (glyco)peptides
- LC-MS analysis
- Data analysis
  - Protein identification (Mascot) and visualization of protein identities ( Scaffold)
  - Filtering for proteins with consensus sequence N-X-S/T (X≠P).
- Expected turnaround time (after arrival of samples):
  - 3-7 days

### Requirements and considerations:

- The enrichment of glycoproteins should be optimized and performed by the customer.
- submit samples as Coomassie Blue stained gel bands (whenever possible avoid silver staining) or in solution (free of protease inhibitors!)

## N-linked glycans site occupancy

### General description:

A glycoprotein may contain a number of potentially N-glycosylated sites, each of which may or may not be glycosylated.

The workflow consists of:

- Glycoprotein enrichment with specific protocols
- N-Glycan release using PNGase F or A (OPTION: in H<sub>2</sub><sup>18</sup>O water)
- Protein digestion (Trypsin) to generate peptides
- LC-MS analysis
- Data analysis
  - Protein identification (Mascot) and visualization of protein identities ( Scaffold)
  - Filtering for consensus sequence N-X-S/T (X≠P) with deamidated N.
- Expected turnaround time (after arrival of samples):
  - 5-10 days

### Requirements and considerations:

- The enrichment of glycoproteins should be optimized and performed by the customer.
- The usage of H<sub>2</sub><sup>18</sup>O water significantly reduces the number of false positive results.
- For the usage of H<sub>2</sub><sup>18</sup>O, please contact us at [proteomics@fgcz.ethz.ch](mailto:proteomics@fgcz.ethz.ch)

## Site-specific glycosylation

### General description:

The N-Glycan microheterogeneity can be studied analysing what kind of glycan is attached to which site of a protein.

The workflow consists of:

- Glycoprotein enrichment with specific protocols
- Protein digestion to generate (glyco)peptides
- LC-MS analysis using different fragmentation techniques
- Data analysis
  - Glycopeptide identification using Mascot or Byonic
  - Manual validation by visual inspection of spectra
- Expected turnaround time (after arrival of samples):
  - 5-15 days

### Requirements and considerations:

- The enrichment of glycoprotein should be optimized and performed by the customer.
- The manual inspection is required, depending on availability of a glycan database for a given species.

## Characterization of N- and O-linked glycans

### General description:

A glycome often comprises a repertoire of closely related structures, many of which are structural isomers. The scope is to know detailed glycan structures or general glycan profiling.

The workflow consists of:

- Glycoprotein enrichment with specific protocols
- N-Glycans: release using PNGase F or A / O-Glycans:  $\beta$ -elimination/reduction reaction
- Chemical derivation, mainly using the permethylation reaction
- MALDI-MS analysis
- Data analysis
  - Glycan identification using ProteinScape 3.0 and Glycoworkbench 2
  - Annotation of the glycan structures and their absolute intensities
- Expected turnaround time (after arrival of samples):
  - 10-15 days

### Requirements and considerations:

- The enrichment of glycoprotein should be optimized and performed by the customer.
- The analysis of N-Glycans can be offered as a service, while the analysis of O-Glycans is currently handled only via User Lab.